

Tips for Optimal Placement of Immobilizer



Harness Position at Opposite Shoulder



Abductor Pillow Placed High in Axilla

Valve Access Panel
Should be positioned
posteriorly with BREG
logo facing forward.

Adjust Harness straps to
suspend pillow and allow
operative arm to rest on
pillow tension free.

Arm and Wrist Velcro
Straps are then
positioned and applied
loosely for support



Pillow Positioned to Keep Elbow Forward of Midline

Once the pillow is placed high in axilla the arm should be draped loosely over pillow for comfort.

Elbow should remain forward of midline when adjusting and applying the arm and wrist holders.



Posterior View

Larger Velcro Loop is used for arm to help prevent it from falling backwards beyond the midline



Side View

Ball may be placed to allow wrist and hand to rest comfortably.

